Návod k obsluze stroje Agilent 6530 Q-TOF LC/MS

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Method Editor window

Spectrum Plot window

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The screenshot shows the Agilent MassHunter Workstation Data Acquisition software interface. The interface is divided into several main sections:

- Instrument Status Window:** Located at the top left, it displays the status of various instrument components: Sampler (Idle), Binary Pump (Standby), Column Comp. (Not Ready), VWD (Not Ready), and Q-TOF (Not Ready). It also shows a status message "Instrument Not Ready" and control buttons for On and Off.
- Chromatogram Plot window:** Located in the middle left, it displays a Total Ion Chromatogram (TIC) plot with a y-axis ranging from -1000000 to 8000000 and an x-axis ranging from 0 to 60 minutes.
- Method Editor window:** Located in the bottom left, it displays the configuration for the current method, including the path, estimated run time, and options for pre-run and post-run scripts.
- Spectrum Plot window:** Located in the middle right, it displays a UV-Vis spectrum plot with a y-axis labeled "mAU" and an x-axis labeled "nm" ranging from -5 to 105.
- Actual window:** Located at the top right, it displays a table of parameters and their values.

| Parameter | Value |
|----------------------------------|---------------|
| Q-TOF: TOF Vac | 1.64E-07 Torr |
| Q-TOF: Quad Vac | 2.43E-05 Torr |
| Q-TOF: Drying Gas | 3.0 l/min |
| Q-TOF: Error State | - |
| Q-TOF: Ready State | False |
| Q-TOF: Rough Vac | 2.81E+00 Torr |
| Q-TOF: Vaporizer/Sheath Gas Temp | 125 °C |
| Q-TOF: Gas Temp | 300 °C |

At the bottom of the interface, there is a navigation bar with three tabs: "Method Editor", "Worklist", and "Sample Run".

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- e) Uveďte **pozici vzorku** v *Sample Position* (2-99).

| | Sample Name | Sample Position | Method | Data File | Sample Type | Level Name | Inj Vol (µl) | Comment | Sample Group | Info. |
|----|--|-----------------|-----------------------------|----------------------------------|-------------|------------|--------------|---------|--------------|-------|
| 1> | blank20.1.2020 | Vial 1 | general_short_positive(400u | D:\MassHunter\Data\Jindrich\Kasa | Sample | | As Method | | | |
| 2 | PK1069 | Vial 61 | general_short_positive(400u | D:\MassHunter\Data\Jindrich\Kasa | Sample | | As Method | | | |
| 3 | washing | Vial 1 | general_short_positive(400u | D:\MassHunter\Data\Jindrich\Kasa | Sample | | As Method | | | |
| 4 | Script: SCP_PumpsAllOff()\{MH_Acq_Scripts.exe} | | | | | | | | | |

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Chromatogram Plot

TIC Binary Pump: Pressure VWD: Signal A

Spectrum Pane

UV | UV

mAU x10⁻²


nm

Worklist

| | Sample Name | Sample Position | Method | Data File | Sample Type | Level Name | Inj Vol (µl) | Comment | Sample Group | Info. |
|---|--|-----------------|-----------------------------|--|-------------|------------|--------------|---------|--------------|-------|
| 1 | blank20.1.2020 | Vial 1 | general_short_positive(400u | D:\MassHunter\Data\Jindrich\Kasal\Sample | | | As Method | | | |
| 2 | PK1069 | Vial 61 | general_short_positive(400u | D:\MassHunter\Data\Jindrich\Kasal\Sample | | | As Method | | | |
| 3 | washing | Vial 1 | general_short_positive(400u | D:\MassHunter\Data\Jindrich\Kasal\Sample | | | As Method | | | |
| 4 | Script: SCP_PumpsAll0#0\MH_Acq_Scripts.exe | | | | | | | | | |

Errors and warnings

1. Skipping overlapped injection because High throughput Optimization is not en...
2. Skipping overlapped injection because High throughput Optimization is not en...
3. Skipping overlapped injection because High throughput Optimization is not en...

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